

REAL-WORLD OUTCOMES WITH WATCHMAN FLX™: EARLY RESULTS FROM SURPASS¹



Late Breaking Clinical Trial, Presented at CRT 2022

Authors: Samir R. Kapadia^A, Robert W Yeh^B, Matthew J Price^C, Jonathan P Piccini^D, Devi G Nair^E, Jonathan C Hsu^F, James V Freeman^G, Thomas Christen^H, Dominic J Allocco^H and Douglas N Gibson^C

A. Department of Cardiovascular Medicine, Heart, Vascular, and Thoracic Institute, Cleveland Clinic, Cleveland, OH B. Beth Israel Deaconess Medical Center, Boston, MA C. Division of Cardiovascular Diseases, Scripps Clinic, La Jolla, CA D. Division of Cardiology, Duke University Medical Center, Durham, NC E. St Bernard's Heart and Vascular Center, Jonesboro, AR F. Cardiac Electrophysiology Section, Division of Cardiology, Department of Medicine, University of California, San Diego, CA G. Section of Cardiovascular Medicine, Yale University School of Medicine, New Haven, Connecticut; Center for Outcomes Research and Evaluation, Yale-New Haven Hospital, New Haven, CT H. Boston Scientific Corp, Marlborough MA

The SURPASS analysis of the NCDR LAAO registry includes the largest number of commercial WATCHMAN FLX patients to date. This data reinforces the outstanding safety of WATCHMAN FLX with a 0.37% major procedural adverse event rate at 7 days or hospital discharge, and a 0.28% ischemic stroke rate through 45 days in > 16,000 real-world NVAF patients.¹

SURPASS Design¹

- SURPASS assesses safety and efficacy outcomes in patients in the NCDR-LAAO Registry™ that received a commercial WATCHMAN FLX device.
- This SURPASS analysis includes 16,048 patients receiving a WATCHMAN FLX device between August 5th, 2020, and March 31st, 2021.
- SURPASS will continue to assess WATCHMAN FLX patients included in the NCDR-LAAO Registry from August 2020 through August 2022 and will follow these patients through 2 years post-implant.

SURPASS Endpoints¹

Safety Endpoint	Composite of all-cause death, ischemic stroke, systemic embolism, or device/procedure-related events requiring open cardiac surgery or major endovascular intervention between device implantation and 7 days or hospital discharge (whichever is later)
Effectiveness Endpoint	Occurrence of ischemic stroke or systemic embolism at 24 months post-implant
Additional Endpoints	<ul style="list-style-type: none"> • All-Cause Death • Stroke • Device-Related Thrombus • Systemic Embolism • Major Bleeding • Effective Device Closure • Implant Success • Device Embolization

Mean baseline patient characteristics: (N=16,048)

- Average Age/years: 76.1 ± 7.9
- CHA₂DS₂-VASc: 4.8 ± 1.5
- HAS-BLED - 2.4 ± 1.0
- 40% women
- 62% with prior clinically relevant bleeding event

SURPASS Outcomes

Key Safety Endpoint within 7 days or Hospital Discharge*¹

Event	WATCHMAN FLX N=16,048		
Key Safety Endpoint	0.37% (60/16,048) 95% CI [0.29, 0.48]		
N=60	Death (n=21)	Ischemic Stroke (n=19)	Events requiring intervention (n=20)

The very low **0.37%** procedural adverse event rate in 16,048 commercial patients confirms the trusted safety profile of WATCHMAN FLX in a broader clinical practice setting.

*Defined as the occurrence of all-cause death, ischemic stroke, systemic embolism, or device or procedure related events requiring open cardiac surgery or major endovascular intervention between device implantation and seven days or hospital discharge (whichever is later)

(SURPASS Outcomes Continued)

SURPASS Outcomes (Continued)


	Discharge N=16,048	45-days N=14,107
Pericardial effusion requiring either surgical or percutaneous intervention	0.32% (52/16,048)	0.51% (72/14,107)
• PE requiring open cardiac surgery	0.01% (2/16,048)	0.03% (4/14,107)
• PE requiring percutaneous treatment	0.31% (50/16,048)	0.50% (70/14,107)

0.32%
of patients had a pericardial effusion requiring intervention at discharge.¹

Excellent procedural success with **98%** of patients implanted (N=16,048/16,446)¹.

These data reinforce the high procedural success and low major procedure-related complication rate of WATCHMAN FLX in this large real-world experience.



Key Clinical Events at 45 days¹



SURPASS data demonstrates low adverse events and a **0.28%** ischemic stroke rate at 45 days.

Comparison with PINNACLE FLX^{2*}

45-day Outcomes



* Results from different clinical investigations are not directly comparable.

These SURPASS Data reinforces the excellent safety profile WATCHMAN FLX demonstrated in the PINNACLE FLX Trial, with largest (n=16,048) real-world WATCHMAN FLX patients studied to date.

References:

¹ Late Breaking Clinical Trial at CRT 2022, Presented by Dr. Samir Kapadia

² Kar, S., Circulation, 2021

Brief Summary

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INTENDED USE/INDICATIONS FOR USE

The WATCHMAN FLX Device is indicated to reduce the risk of thromboembolism from the left atrial appendage in patients with non-valvular atrial fibrillation who:

- Are at increased risk for stroke and systemic embolism based on CHADS₂ or CHA₂DS₂-VASc scores and are recommended for anticoagulation therapy;
- Are deemed by their physicians to be suitable for anticoagulation therapy; and
- Have an appropriate rationale to seek a non-pharmacologic alternative to anticoagulation therapy, taking into account the safety and effectiveness of the device compared to anticoagulation therapy.

CONTRAINDICATIONS

Do not use the WATCHMAN FLX Device if:

- Intracardiac thrombus is present.
- An atrial septal defect repair or closure device or a patent foramen ovale repair or closure device is present.
- The LAA anatomy will not accommodate a Closure Device (see **Table 45 of the eIFU**).
- The patient has a known hypersensitivity to any portion of the device material or the individual components (see Device Description section of the eIFU) such that the use of the WATCHMAN FLX Device is contraindicated.
- Any of the customary contraindications for other percutaneous catheterization procedure (e.g., patient size too small to accommodate TEE probe or required catheters) or conditions (e.g., active infection, bleeding disorder) are present.
- There are contraindications to the use of anticoagulation therapy, aspirin, or P2Y₁₂ inhibitor.

WARNINGS

Implantation of the WATCHMAN FLX Device should only be performed by interventional cardiologists and/ or electrophysiologists who are trained in percutaneous and transseptal procedures and who have completed the WATCHMAN FLX Physician Training program.

- This device has not been studied in pregnant or breastfeeding women. Careful consideration should be given to use of the Closure Device in pregnant and/ or breastfeeding women due to the risk of significant exposure to x-rays and the use of anticoagulation medication.
- Device selection should be based on accurate LAA measurements obtained using echocardiographic imaging guidance in multiple views (TEE recommended in multiple angles [e.g., 0°, 45°, 90°, 135°]) to avoid improper Closure Device sizing.
- Do not release (i.e., unscrew) the WATCHMAN FLX Device from the core wire unless all release criteria are satisfied to avoid suboptimal results.
- Potential for Closure Device embolization exists with cardioversion < 30 days following Closure Device implantation; verify Closure Device position after cardioversion during this period.
- Appropriate post-procedure drug therapy should be followed. See Post-Procedure Information section (of the eIFU) for further detail.

PRECAUTIONS

- The safety and effectiveness (and benefit-risk profile) of the WATCHMAN FLX Device has not been established in patients for whom long-term anticoagulation is determined to be contraindicated.
- The LAA is a thin-walled structure. Use caution when accessing the LAA, and deploying, recapturing, and repositioning the Closure Device.
- Use caution when introducing a WATCHMAN Access System to prevent damage to cardiac structures.
- Use caution when introducing the Delivery System to prevent damage to cardiac structures.

- To prevent damage to the Delivery Catheter or Closure Device, do not allow the WATCHMAN FLX Device to protrude beyond the distal tip of the Delivery Catheter when inserting the Delivery System into the Access Sheath.
- If using a power injector, the maximum pressure should not exceed 100 psi.

PATIENT SELECTION FOR TREATMENT

In considering the use of the WATCHMAN FLX Device, the rationale for seeking an alternative to long-term anticoagulation therapy and the safety and effectiveness of the device compared to anticoagulation should be taken into account.

- The presence of indication(s) for long-term anticoagulation therapy, other than non-valvular atrial fibrillation (e.g. mechanical heart valve, hypercoagulable states, recurrent deep venous thrombosis).

Details regarding the indications, contraindications, warnings, and precautions for oral anticoagulants approved for patients with non-valvular atrial fibrillation are provided in their respective Instructions for Use.

Of note:

- The safety and effectiveness (and benefit-risk profile) of the WATCHMAN FLX Device has not been established in patients for whom long-term anticoagulation is determined to be contraindicated.

Factors that need to be considered for the WATCHMAN FLX Device and implantation procedure include the following:

- Overall medical status, including conditions which might preclude the safety of a percutaneous, transcatheter procedure.
- Suitability for percutaneous, transseptal procedures, including considerations of:
 - Cardiac anatomy relating to the LAA size and shape.
 - Vascular access anatomy (e.g., femoral vein size, thrombus, or tortuosity).
 - Ability of the patient to tolerate general or local anesthesia.
 - Ability of the patient to undergo required imaging.
- Ability to comply with the recommended post-WATCHMAN FLX Device implant pharmacologic regimen (see Post-Procedure Information section) especially for patients at high risk for bleeding.

ADVERSE EVENTS

Potential adverse events (in alphabetical order) which may be associated with the use of a left atrial appendage closure device or implantation procedure include but are not limited to:

Air embolism, Airway trauma, Allergic reaction to the contrast media, anesthetic, WATCHMAN Implant material, or medications, Altered mental status, Anemia requiring transfusion, Anesthesia risks, Angina, Anoxic encephalopathy, Arrhythmias, Atrial septal defect, Bruising, hematoma, or seroma near the catheter insertion site, Cardiac perforation, Chest pain/ discomfort, Confusion post procedure, Congestive heart failure, Contrast related nephropathy, Cranial bleed, Death, Decreased hemoglobin, Deep vein thrombosis, Device embolism, Device fracture, Device thrombosis, Edema, Embolism, Excessive bleeding, Fever, Fistula, Groin pain, Groin puncture bleed, Hematuria, Hemoptysis, Hypotension, Hypoxia, Improper wound healing, Inability to reposition, recapture, or retrieve the device, Infection/pneumonia, Interatrial septum thrombus, Intratracheal bleeding, Major bleeding requiring transfusion, Misplacement of the device/ improper seal of the appendage/movement of device from appendage wall, Myocardial erosion, Nausea, Oral bleeding, Pericardial effusion/ tamponade, Pleural effusion, Prolonged bleeding from a laceration, Pulmonary edema, Renal failure, Respiratory insufficiency/failure, Stroke - Hemorrhagic, Stroke - Ischemic, Surgical removal of the device, TEE complications (e.g., throat pain, bleeding, esophageal trauma), Thrombocytopenia, Thrombosis, Transient ischemic attack (TIA), Valvular or vascular damage, Vasovagal reactions

There may be other potential adverse events that are unforeseen at this time.

92574167 A.1

CE 2797

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.