


FAST TRACK

to your therapy delivery
in a **SINGLE SOLUTION**




VersaCross Connect™

LAAC Access Solution

DESIGNED for use with
WATCHMAN
Access Systems†


**3-IN-1
RF WIRE**
for ZERO EXCHANGE delivery
of **WATCHMAN™** Access Sheath†



Smooth transition between
VersaCross Connect™ Dilator and
WATCHMAN™ Access Sheath†



Snap-lock designed for
WATCHMAN™ Access Sheaths†



*The therapy referred to here may be one of several left heart procedures that require transseptal access. Additional procedural steps and/or devices may be required to deliver the therapy. Before use, consult Instructions for Use for any devices accordingly.

† WATCHMAN™ Access Sheaths as shown in this brochure is sold separately.

‡ The VersaCross Connect™ Transseptal Dilator is for use with a 12F (4.2 mm) ID WATCHMAN™ Access Sheath that is 75 cm in length or the WATCHMAN TruSteer Access Sheath that is 67 cm in length, specifically: WATCHMAN™ Access System [Models: M635TU40060, M635TU10060, M635TU20060]; WATCHMAN™ TruSeal™ Access System [Models: M635TU70010, M635TU70040, M635TU70020]; WATCHMAN FXD Curve™ Access System [Models: M635TU80010, M635TU80020]; WATCHMAN TruSteer™ Access System [Model: M635TU90050].

All trademarks are property of their respective owners. Patents Pending and/or issued. Caution: Federal law (USA) restricts these devices to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse events, and Operator's instructions. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries.



For brief summaries,
scan QR code or visit
qrco.de/bs-vxawm

**Boston
Scientific**
Advancing science for life™

Baylis Medical Company Inc.
5959 Trans-Canada Highway
Montreal, QC Canada H4T 1A1
www.bostonscientific.com

Customer Service
United States: 1-888-272-1001
Canada: 1-888-359-9691

© 2025 Boston Scientific Corporation
or its affiliates. All rights reserved.
EP-1536805-AC

**Boston
Scientific**
Advancing science for life™


VersaCross Connect™
LAAC Access Solution

ZERO EXCHANGE
WATCHMAN™ Access Systems

ZERO EXCHANGE

WATCHMAN™ Access Sheath Delivery⁺

EXCHANGELESS 3-in-1 RF Wire



1 Start at the SVC



2 RF transseptal puncture



3 Deliver WATCHMAN™ Sheath

INSTANTLY DEPLOY


0.035" WIRE

to deliver WATCHMAN™ Access Systems⁺ with confidence


DELIVER THERAPY ON TARGET



TRUform™
Shapeable Technology




Precision RF
Puncture Technology



OMNIvizi™ Technology

Reliably locate the VersaCross Connect solution on fluoroscopy and ultrasound. Positional markers indicate position of RF tip within dilator.

the right **CURVE** at the right **TIME**



TRANSSEPTAL CURVE

VersaCross Connect™ Transseptal Dilator with TRUform™ shapeable technology optimizes transseptal approach



THERAPY DELIVERY CURVE

WATCHMAN™ Access Sheaths retain their curve for therapy delivery

Personalize your solution

1. Choose your VersaCross™ RF Wire

2. Choose your VersaCross Connect™ Transseptal Dilator



Length: 180 cm, 230 cm
Diameter: 0.035"



Inner Diameter: 0.035"
Outer Diameter: 12F (4.2 mm)
Useable Length: 85 cm

Solution with Pigtail Wire

Part No.	Dilator curve	Wire length (cm)
VXAK0101	D0	180
VXAK0103	D1	180
VXAK0105	D0	230
VXAK0107	D1	230

Solution with J-Tip Wire

Part No.	Dilator curve	Wire length (cm)
VXAK0109	D0	180