



EASY. PRECISE. OPTIMIZED FOR WATCHMAN FLX.™

WATCHMAN TruPlan is a CT imaging software that enhances patient screening, procedure planning, and accurate diagnosis for LAAC.

DEVELOPED FOR YOU BY MARKET LEADERS BOSTON SCIENTIFIC AND CIRCLE CVI

Minimum Specifications (v3.1)

Below are system requirements for standalone, desktop client and server installations.

Requirement	Standalone Installation	Desktop Client	Server
CPU	Windows: Intel i5 7th Gen Mac: Intel Core i5 7th gen+, Apple silicon	Windows: Intel i5 7th Gen Mac: Intel Core i5 7th gen+, Apple silicon	4-core Xeon
RAM	16GB	8GB	16GB
VIDEO	Intel Integrated HD Graphics 630+	Intel Integrated HD Graphics 630+	n/a
MAC	macOS 11, 12	macOS 11, 12	n/a
WINDOWS	Windows 10, 11	Window 10, 11	Windows Server 2016 or newer
Display	1920x1080	1920x1080	n/a

**It is advised to download and install the latest graphics driver when using TruPlan.*

EASY. PRECISE. OPTIMIZED FOR WATCHMAN FLX.

Platform Functionality

- **Technology Platform** – HIPAA-compliant software with standalone or client/server deployment
- **System Access** – Software-specific login credentials are not required; Software-specific login credentials are not required; user management is enabled via LDAP integration
- **Data Used** – DICOM data containing PHI is imported directly from PACS or locally (via hard drive or storage device)
- **Data Access** – Data is secured within existing IT infrastructure; the software does not transmit PHI to the vendor
- **Data Connectivity** – Optional connectivity to PACS using standard DICOM protocols
- **Internet Connectivity** – The software runs locally without internet connection; license key activation will require a direct internet connection

Cybersecurity



PHI

To protect patient information, it is advised to make sure the operating system of the computer is protected with a password and the patient information is stored on an encrypted disk partition.



SECURITY

WATCHMAN TruPlan is tested with the security suite of Windows Defender, there are no known limitations. It is advised to install the latest security updates for the operating system.



BACKUPS

WATCHMAN TruPlan is not meant as the single storage location of the DICOM data, therefore there are no requirements for making backups of the system.

NEED HELP? CONTACT CIRCLE CVI TECH SUPPORT

CALL +1 (403) 338-1870

VISIT www.watchman.com/truplan

EMAIL support@circlecvi.com