



BUILT FOR **CONFIDENCE**. BUILT FOR **COLLABORATION**. BUILT FOR **LAAC**.

WATCHMAN TruPlan is a CT imaging software to facilitate patient screening and procedure planning for LAAC


DEVELOPED FOR YOU BY MARKET LEADERS BOSTON SCIENTIFIC AND CIRCLE CVI

Minimum Specifications (v2.0)




OPERATING SYSTEM

Windows 10 64bit, Mac OS 10.15 or 11.x



PROCESSOR

Intel Core i5 or better



SYSTEM RAM

8 GB (16 GB recommended)



VIDEO

Intel HD 620 or better (Dedicated GPU recommended)*



DISPLAY

1280 x 800

**It is advised to download and install the latest graphics driver when using TruPlan.*

BUILT FOR **CONFIDENCE**. BUILT FOR **COLLABORATION**. BUILT FOR **LAAC**.

Platform Functionality

- **Technology Platform** – HIPAA-compliant standalone software; not server or cloud based
- **System Access** – Software-specific login credentials are not required; system access is controlled by the workstation login credentials where the software is installed
- **Data Used** – Dicom data containing PHI is imported directly from PACS or locally (via hard drive or storage device)
- **Data Access** – Data is secured within existing IT infrastructure; the software does not transmit PHI to the vendor
- **Data Connectivity** – Optional connectivity to PACS using Orthanc, an open-source, lightweight DICOM server
- **Internet Connectivity** – The software runs locally without internet connection; license key activation will require a direct internet connection

Cybersecurity



PHI

To protect patient information, it is advised to make sure the operating system of the computer is protected with a password and the patient information is stored on an encrypted disk partition.



SECURITY

TruPlan is tested with the security suite of Windows Defender, there are no known limitations. It is advised to install the latest security updates for the operating system.



BACKUPS

TruPlan is not meant as the single storage location of the DICOM data, therefore there are no requirements for making backups of the system.

NEED HELP? CONTACT CIRCLE CVI TECH SUPPORT

CALL +1 (403) 338-1870

VISIT www.watchman.com/truplan

EMAIL support@circlecvi.com

Brief Summary: Indications, contraindications, warnings and precautions can be found in the product labeling. TruPlan shall be used by qualified medical professionals, experienced in examining and evaluating medical images.

©2022 Boston Scientific Corporation or its affiliates. All rights reserved. All trademarks are property of their respective owners. WATCHMAN TruPlan software is developed and owned by Circle Cardiovascular Imaging Inc. (Calgary, AB, Canada), and Boston Scientific is the exclusive reseller of WATCHMAN TruPlan software. SH-866612-AD

Please visit www.circlecvi.com/watchman-truplan for a summary of safety information and to download the full Instructions for Use (User Manual). Rx Only.